

von Willebrand Factor

Bio-Immuno-Assays for determination of von Willebrand Factor Collagen Binding Activity and von Willebrand Factor Antigen

technoclone

A clear view into the future!

TECHNOZYM® vWF:CBA AND TECHNOZYM® vWF:AG

The von Willebrand disease (vWD) is a bleeding disorder resulting from a qualitative or quantitative deficiency of the von Willebrand factor (vWF).

To detect and distinguish between quantitative vWF deficiencies in types 1 and 3 vWD and qualitative vWF deficiencies in types 2A, 2B and 2M vWD, the collagen binding assay (vWF:CBA), which measures the binding of preferentially High Molecular Weight vWF multimers to collagen, is combined with vWF antigen determination.

Full pannel of assays for accurate diagnosis of vWD

TECHNOZYM[®] vWF:Ag

ELISA for quantitative determination of vWF antigen

TECHNOZYM® vWF:CBA

ELISA coated with human Collagen Type III for detection of qualitative vWF deficiencies in types 2A and 2B vWD

◆ TECHNOZYM[®] vWF:CBA Collagen Type I

ELISA coated with human Collagen Type I for detection of qualitative vWF deficiencies in types 2A and 2B vWD

TECHNOZYM[®] vWF:CBA Collagen Type VI

ELISA coated with human Collagen Type IV is recommended for an accurate diagnosis of type 2M vWD

Product	Contents	REF	Package
TECHNOZYM® vWF:Ag ELISA Kit	12 x 8 anti-vWF strips 1 x 0.3 ml anti-vWF Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450201	96 tests
TECHNOZYM® vWF:CBA ELISA Kit	12 x 8 Collagen Type III strips 1 x 0.3 ml anti-vWF Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450301	96 tests
TECHNOZYM [®] vWF:CBA Collagen Type I ELISA Kit	12 x 8 Collagen Type I strips 1 x 0.3 ml anti-vWF Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450311	96 tests
TECHNOZYM® vWF:CBA Collagen Typ VI ELISA Kit	12 x 8 Collagen Type VI strips 1 x 0.3 ml anti-vWF Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450321	96 tests

Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH

Brunner Str. 67 | 1230 Vienna | Austria Phone: +43 1 86373-0 | Fax: +43 1 86373-44 products@technoclone.com | www.technoclone.com













